# IMPRELON® »S« – Hinweise / Hints

0,5 mm - 0,75 mm - 1,0 mm - 1,5 mm - 2,0 mm - 3,0 mm

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IMPRELON® »S« in the a. m. thicknesses is a high quality polycarbonate which is a break- and abrasionresistant clear thermoplastic. The material is best suited for all splint therapies.

**Pre-drying:** Due to the properties of polycarbonate and the high pressurizing temperature of approximately 180 °C **pre-drying** of IMPRELON® »S« is absolutely **necessary** to **avoid air bubbles** when plastifying this material. IMPRELON® »S« can be stored in the BIOSTAR-drying oven or, depending on the thickness of the material can be **pre-dried** for the recommended conditioning time in the BIOSTAR drying oven or a similar oven at a temperature of **80 °C – 100 °C**. It may be necessary to prolong first pre-drying time of 2.0 and 3.0 mm sheets depending on humidity, storing and transport conditions. Sheets prepared in this way can be used for approximately 7 days without any air bubbles. Repeated pre-drying is possible without any loss of quality.

**Bite splints:** The etching liquid PERLIBOND can be applied to IMPRELON® »S« if resin should be added to bite splint. For this use 1 – 2 times PERLIBOND for the occlusal areas and model the desired shape by applying autopolymerisate (e. g. DURASPLINT®, STEADY-RESIN) on the moistened bonding agent.

Polymerization and polishing is done as usual.

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IMPRELON® »S« – as well as all our pressure moulding foils – has been tested at the University of Heidelberg with regard to toxic properties, gingiva irritations and mutagenity. The results clearly state that these foils are biologically well-accepted.

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For cleaning of splints made of IMPRELON® »S« we recommend our CETRON®-cleansing and care products.

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